

Calibrating Unit

New Standard Solution Report

Equipment Alcotest 7110 MKIII-C Serial No.: ARTL-0005
Location: PLAINSBORO TOWNSHIP PD
Calibration File No.: 00706 **Calib. Date:** 12/10/2009 **Calib. No.:** 00013
Certification File No.: 00707 **Cert. Date:** 12/10/2009 **Cert. No.:** 00009
Linearity File No.: 00708 **Lin. Date:** 12/10/2009 **Lin. No.:** 00009
Solution File No.: 00724 **Soln. Date:** 01/05/2010 **Soln. No.:** 00095
Sequential File No.: 00724 **File Date:** 01/05/2010

Calibrating Unit: WET **Model No.:** CU-34 **Serial No.:** DDUF S3-0065
Control Solution %: 0.100% **Expires:** 07/23/2010
Solution Control Lot: 08G058 **Bottle No.:** 0921

Function	Result %BAC	Time HH:MM	Temperature Simulator (°C)	Comment(s) or Error(s)
Ambient Air Blank	0.000%	12:30S		
Control 1 EC	0.100%	12:31S	34.0°C	*** TEST PASSED ***
Control 1 IR	0.101%	12:31S	34.0°C	*** TEST PASSED ***
Ambient Air Blank	0.000%	12:32S		
Control 2 EC	0.099%	12:32S	34.0°C	*** TEST PASSED ***
Control 2 IR	0.099%	12:32S	34.0°C	*** TEST PASSED ***
Ambient Air Blank	0.000%	12:33S		
Control 3 EC	0.099%	12:34S	34.1°C	*** TEST PASSED ***
Control 3 IR	0.099%	12:34S	34.1°C	*** TEST PASSED ***
Ambient Air Blank	0.000%	12:34S		

All tests within acceptable tolerance.

On this date, I installed the above indicated "NEW SOLUTION" in accordance with Alcotest 7110 operator training and procedures established by the (NJSP) Chief Forensic Scientist.

Changed By:

Last Name: JANKOWSKI

First Name: JOSEPH

MI: E

Signature: 

Badge No.: 714

Date: 01/05/2010

Calibrating Unit

New Standard Solution Report

Equipment Alcotest 7110 MKIII-C Serial No.: ARTL-0005
Location: PLAINSBORO TOWNSHIP PD
Calibration File No.: 00706 **Calib. Date:** 12/10/2009 **Calib. No.:** 00013
Certification File No.: 00707 **Cert. Date:** 12/10/2009 **Cert. No.:** 00009
Linearity File No.: 00708 **Lin. Date:** 12/10/2009 **Lin. No.:** 00009
Solution File No.: 00740 **Soln. Date:** 01/25/2010 **Soln. No.:** 00096
Sequential File No.: 00740 **File Date:** 01/25/2010

Calibrating Unit: WET **Model No.:** CU-34 **Serial No.:** DDUF S3-0065
Control Solution %: 0.100% **Expires:** 07/23/2010
Solution Control Lot: 08G058 **Bottle No.:** 0923

Function	Result	Time	Temperature	Comment(s)
	%BAC	HH:MM	Simulator (°C)	or Error(s)
Ambient Air Blank	0.000%	08:33S		
Control 1 EC	0.099%	08:34S	34.1°C	*** TEST PASSED ***
Control 1 IR	0.100%	08:34S	34.1°C	*** TEST PASSED ***
Ambient Air Blank	0.000%	08:35S		
Control 2 EC	0.098%	08:35S	34.1°C	*** TEST PASSED ***
Control 2 IR	0.100%	08:35S	34.1°C	*** TEST PASSED ***
Ambient Air Blank	0.000%	08:36S		
Control 3 EC	0.099%	08:37S	34.1°C	*** TEST PASSED ***
Control 3 IR	0.100%	08:37S	34.1°C	*** TEST PASSED ***
Ambient Air Blank	0.000%	08:37S		

All tests within acceptable tolerance.


On this date, I installed the above indicated "NEW SOLUTION" in accordance with Alcotest 7110 operator training and procedures established by the (NJSP) Chief Forensic Scientist.

Changed By:

Last Name: JANKOWSKI

First Name: JOSEPH

MI: E

Signature: 

Badge No.: 714

Date: 01/25/2010

Calibrating Unit

New Standard Solution Report

Equipment	Alcotest 7110 MKIII-C	Serial No.:	ARTL-0005
Location:	PLAINSBORO TOWNSHIP PD		
Calibration File No.:	00706	Calib. Date:	12/10/2009
Certification File No.:	00707	Cert. Date:	12/10/2009
Linearity File No.:	00708	Lin. Date:	12/10/2009
Solution File No.:	00746	Soln. Date:	02/24/2010
Sequential File No.:	00746	File Date:	02/24/2010
Calibrating Unit:	WET	Model No.:	CU-34
Control Solution %:	0.100%	Serial No.:	DDUF S3-0065
Solution Control Lot:	08G058	Expires:	07/23/2010
		Bottle No.:	0924

Function	Result %BAC	Time HH:MM	Temperature Simulator (°C)	Comment(s) or Error(s)
Ambient Air Blank	0.000%	18:02S		
Control 1 EC	0.099%	18:02S	34.1°C	*** TEST PASSED ***
Control 1 IR	0.099%	18:02S	34.1°C	*** TEST PASSED ***
Ambient Air Blank	0.000%	18:03S		
Control 2 EC	0.098%	18:04S	34.1°C	*** TEST PASSED ***
Control 2 IR	0.099%	18:04S	34.1°C	*** TEST PASSED ***
Ambient Air Blank	0.000%	18:04S		
Control 3 EC	0.099%	18:05S	34.1°C	*** TEST PASSED ***
Control 3 IR	0.099%	18:05S	34.1°C	*** TEST PASSED ***
Ambient Air Blank	0.000%	18:06S		

All tests within acceptable tolerance.


On this date, I installed the above indicated "NEW SOLUTION" in accordance with Alcotest 7110 operator training and procedures established by the (NJSP) Chief Forensic Scientist.

Changed By:

Last Name: JANKOWSKI

First Name: JOSEPH

MI: E

Signature: 

Badge No.: 714

Date: 02/24/2010

Calibrating Unit

New Standard Solution Report

Equipment Alcotest 7110 MKIII-C Serial No.: ARTL-0005
Location: PLAINSBORO TOWNSHIP PD
Calibration File No.: 00706 Calib. Date: 12/10/2009 Calib. No.: 00013
Certification File No.: 00707 Cert. Date: 12/10/2009 Cert. No.: 00009
Linearity File No.: 00708 Lin. Date: 12/10/2009 Lin. No.: 00009
Solution File No.: 00754 Soln. Date: 03/17/2010 Soln. No.: 00098
Sequential File No.: 00754 File Date: 03/17/2010

Calibrating Unit: WET Model No.: CU-34 Serial No.: DDUF S3-0065
Control Solution %: 0.100% Expires: 07/23/2010
Solution Control Lot: 08G058 Bottle No.: 0925

Function	Result %BAC	Time HH:MM	Temperature Simulator (°C)	Comment(s) or Error(s)
Ambient Air Blank	0.000%	07:43S		
Control 1 EC	0.100%	07:44S	34.1°C	*** TEST PASSED ***
Control 1 IR	0.099%	07:44S	34.1°C	*** TEST PASSED ***
Ambient Air Blank	0.000%	07:45S		
Control 2 EC	0.098%	07:45S	34.0°C	*** TEST PASSED ***
Control 2 IR	0.100%	07:45S	34.0°C	*** TEST PASSED ***
Ambient Air Blank	0.000%	07:46S		
Control 3 EC	0.099%	07:47S	34.0°C	*** TEST PASSED ***
Control 3 IR	0.099%	07:47S	34.0°C	*** TEST PASSED ***
Ambient Air Blank	0.000%	07:47S		

All tests within acceptable tolerance.

On this date, I installed the above indicated "NEW SOLUTION" in accordance with Alcotest 7110 operator training and procedures established by the (NJSP) Chief Forensic Scientist.

Changed By:

Last Name: TAVENER

First Name: FREDERICK

MI: W

Signature: 

Badge No.: 730

Date: 03/17/2010

Calibrating Unit

New Standard Solution Report

Equipment	Alcotest 7110 MKIII-C	Serial No.:	ARTL-0005
Location:	PLAINSBORO TOWNSHIP PD		
Calibration File No.:	00706	Calib. Date:	12/10/2009
Certification File No.:	00707	Cert. Date:	12/10/2009
Linearity File No.:	00708	Lin. Date:	12/10/2009
Solution File No.:	00757	Soln. Date:	03/28/2010
Sequential File No.:	00757	File Date:	03/28/2010
Calibrating Unit:	WET	Model No.:	CU-34
Control Solution %:	0.100%	Serial No.:	DDUF S3-0065
Solution Control Lot:	09L072	Expires:	12/31/2011
		Bottle No.:	0251

Function	Result %BAC	Time HH:MM	Temperature Simulator (°C)	Comment(s) or Error(s)
Ambient Air Blank	0.000%	05:39S		
Control 1 EC	0.100%	05:40S	34.0°C	*** TEST PASSED ***
Control 1 IR	0.100%	05:40S	34.0°C	*** TEST PASSED ***
Ambient Air Blank	0.000%	05:40S		
Control 2 EC	0.100%	05:41S	34.1°C	*** TEST PASSED ***
Control 2 IR	0.100%	05:41S	34.1°C	*** TEST PASSED ***
Ambient Air Blank	0.000%	05:42S		
Control 3 EC	0.100%	05:42S	34.0°C	*** TEST PASSED ***
Control 3 IR	0.101%	05:42S	34.0°C	*** TEST PASSED ***
Ambient Air Blank	0.000%	05:43S		

All tests within acceptable tolerance.

On this date, I installed the above indicated "NEW SOLUTION" in accordance with Alcotest 7110 operator training and procedures established by the (NJSP) Chief Forensic Scientist.

Changed By:

Last Name: JANKOWSKI

First Name: JOSEPH

MI: E

Signature:

 #714

Badge No.: 714

Date: 03/28/2010

Calibrating Unit

New Standard Solution Report

Equipment Alcotest 7110 MKIII-C Serial No.: ARTL-0005
Location: PLAINSBORO TOWNSHIP PD
Calibration File No.: 00706 **Calib. Date:** 12/10/2009 **Calib. No.:** 00013
Certification File No.: 00707 **Cert. Date:** 12/10/2009 **Cert. No.:** 00009
Linearity File No.: 00708 **Lin. Date:** 12/10/2009 **Lin. No.:** 00009
Solution File No.: 00764 **Soln. Date:** 04/23/2010 **Soln. No.:** 00100
Sequential File No.: 00764 **File Date:** 04/23/2010

Calibrating Unit: WET **Model No.:** CU-34 **Serial No.:** DDUF S3-0065
Control Solution %: 0.100% **Expires:** 12/31/2011
Solution Control Lot: 09L072 **Bottle No.:** 0252

Function	Result	Time	Temperature	Comment(s)
	%BAC	HH:MM	Simulator (°C)	or Error(s)
Ambient Air Blank	0.000%	16:44D		
Control 1 EC	0.100%	16:45D	34.1°C	*** TEST PASSED ***
Control 1 IR	0.100%	16:45D	34.1°C	*** TEST PASSED ***
Ambient Air Blank	0.000%	16:46D		
Control 2 EC	0.099%	16:46D	34.1°C	*** TEST PASSED ***
Control 2 IR	0.100%	16:46D	34.1°C	*** TEST PASSED ***
Ambient Air Blank	0.000%	16:47D		
Control 3 EC	0.099%	16:47D	34.1°C	*** TEST PASSED ***
Control 3 IR	0.102%	16:47D	34.1°C	*** TEST PASSED ***
Ambient Air Blank	0.000%	16:48D		

All tests within acceptable tolerance.

On this date, I installed the above indicated "NEW SOLUTION" in accordance with Alcotest 7110 operator training and procedures established by the (NJSP) Chief Forensic Scientist.

Changed By:

Last Name: TAVENER

First Name: FREDERICK

MI: W

Signature: SGT F TAVENER #730

Badge No.: 730

Date: 04/23/2010

Calibrating Unit

New Standard Solution Report

Equipment Alcotest 7110 MKIII-C Serial No.: ARTL-0005
Location: PLAINSBORO TOWNSHIP PD
Calibration File No.: 00706 Calib. Date: 12/10/2009 Calib. No.: 00013
Certification File No.: 00707 Cert. Date: 12/10/2009 Cert. No.: 00009
Linearity File No.: 00708 Lin. Date: 12/10/2009 Lin. No.: 00009
Solution File No.: 00770 Soln. Date: 05/14/2010 Soln. No.: 00101
Sequential File No.: 00770 File Date: 05/14/2010

Calibrating Unit: WET Model No.: CU-34 Serial No.: DDUF S3-0065
Control Solution %: 0.100% Expires: 12/31/2011
Solution Control Lot: 09L072 Bottle No.: 0253

Function	Result %BAC	Time HH:MM	Temperature Simulator (°C)	Comment(s) or Error(s)
Ambient Air Blank	0.000%	15:07D		
Control 1 EC	0.100%	15:07D	34.1°C	*** TEST PASSED ***
Control 1 IR	0.100%	15:07D	34.1°C	*** TEST PASSED ***
Ambient Air Blank	0.000%	15:08D		
Control 2 EC	0.099%	15:08D	34.1°C	*** TEST PASSED ***
Control 2 IR	0.100%	15:08D	34.1°C	*** TEST PASSED ***
Ambient Air Blank	0.000%	15:09D		
Control 3 EC	0.099%	15:10D	34.1°C	*** TEST PASSED ***
Control 3 IR	0.099%	15:10D	34.1°C	*** TEST PASSED ***
Ambient Air Blank	0.000%	15:10D		

All tests within acceptable tolerance.

On this date, I installed the above indicated "NEW SOLUTION" in accordance with Alcotest 7110 operator training and procedures established by the (NJSP) Chief Forensic Scientist.


Changed By:

Last Name: TAVENER

First Name: FREDERICK

MI: W

Signature: _____



Badge No.: 730

Date: 05/14/2010

Calibrating Unit

New Standard Solution Report

Equipment Alcotest 7110 MKIII-C Serial No.: ARTL-0005
Location: PLAINSBORO TOWNSHIP PD
Calibration File No.: 00780 Calib. Date: 06/07/2010 Calib. No.: 00014
Certification File No.: 00781 Cert. Date: 06/07/2010 Cert. No.: 00010
Linearity File No.: 00782 Lin. Date: 06/07/2010 Lin. No.: 00010
Solution File No.: 00783 Soln. Date: 06/07/2010 Soln. No.: 00102
Sequential File No.: 00783 File Date: 06/07/2010

Calibrating Unit: WET Model No.: CU-34 Serial No.: DDUF S3-0065
Control Solution %: 0.100% Expires: 12/31/2011
Solution Control Lot: 09L072 Bottle No.: 0254

Function	Result %BAC	Time HH:MM	Temperature Simulator (°C)	Comment(s) or Error(s)
Ambient Air Blank	0.000%	18:51D		
Control 1 EC	0.101%	18:52D	34.0°C	*** TEST PASSED ***
Control 1 IR	0.101%	18:52D	34.0°C	*** TEST PASSED ***
Ambient Air Blank	0.000%	18:52D		
Control 2 EC	0.101%	18:53D	34.0°C	*** TEST PASSED ***
Control 2 IR	0.101%	18:53D	34.0°C	*** TEST PASSED ***
Ambient Air Blank	0.000%	18:54D		
Control 3 EC	0.101%	18:54D	34.0°C	*** TEST PASSED ***
Control 3 IR	0.102%	18:54D	34.0°C	*** TEST PASSED ***
Ambient Air Blank	0.000%	18:55D		

All tests within acceptable tolerance.

On this date, I installed the above indicated "NEW SOLUTION" in accordance with Alcotest 7110 operator training and procedures established by the (NJSP) Chief Forensic Scientist.

TEMPERATURE PROBE SERIAL # DDUF P2-144 (M)

Changed By:

Last Name: DENNIS

First Name: MARC

MI: W.

Signature: TPL J [Signature] #5925

Badge No.: 5925

Date: 06/07/2010

Calibrating Unit

New Standard Solution Report

Equipment Alcotest 7110 MKIII-C Serial No.: ARTL-0005
Location: PLAINSBORO TOWNSHIP PD
Calibration File No.: 00780 **Calib. Date:** 06/07/2010 **Calib. No.:** 00014
Certification File No.: 00781 **Cert. Date:** 06/07/2010 **Cert. No.:** 00010
Linearity File No.: 00782 **Lin. Date:** 06/07/2010 **Lin. No.:** 00010
Solution File No.: 00792 **Soln. Date:** 06/21/2010 **Soln. No.:** 00103
Sequential File No.: 00792 **File Date:** 06/21/2010

Calibrating Unit: WET **Model No.:** CU-34 **Serial No.:** DDUF S3-0064
Control Solution %: 0.100% **Expires:** 12/31/2011
Solution Control Lot: 09L072 **Bottle No.:** 0260

Function	Result %BAC	Time HH:MM	Temperature Simulator (°C)	Comment(s) or Error(s)
Ambient Air Blank	0.000%	16:11D		
Control 1 EC	0.101%	16:12D	34.0°C	*** TEST PASSED ***
Control 1 IR	0.101%	16:12D	34.0°C	*** TEST PASSED ***
Ambient Air Blank	0.000%	16:12D		
Control 2 EC	0.100%	16:13D	34.0°C	*** TEST PASSED ***
Control 2 IR	0.102%	16:13D	34.0°C	*** TEST PASSED ***
Ambient Air Blank	0.000%	16:13D		
Control 3 EC	0.101%	16:14D	34.0°C	*** TEST PASSED ***
Control 3 IR	0.100%	16:14D	34.0°C	*** TEST PASSED ***
Ambient Air Blank	0.000%	16:15D		

All tests within acceptable tolerance.

On this date, I installed the above indicated "NEW SOLUTION" in accordance with Alcotest 7110 operator training and procedures established by the (NJSP) Chief Forensic Scientist.


TEMPERATURE PROBE SERIAL # DDSF P2-0396

Changed By:

Last Name: DENNIS

First Name: MARC

MI: W.

Signature: 

#5925

Badge No.: 5925

Date: 06/21/2010

Calibrating Unit

New Standard Solution Report

Equipment	Alcotest 7110 MKIII-C	Serial No.:	ARTL-0005
Location:	PLAINSBORO TOWNSHIP PD		
Calibration File No.:	00780	Calib. Date:	06/07/2010
Certification File No.:	00781	Cert. Date:	06/07/2010
Linearity File No.:	00782	Lin. Date:	06/07/2010
Solution File No.:	00808	Soln. Date:	07/17/2010
Sequential File No.:	00808	File Date:	07/17/2010
Calibrating Unit:	WET	Model No.:	CU-34
Control Solution %:	0.100%	Serial No.:	DDUF S3-0064
Solution Control Lot:	09L072	Expires:	12/31/2011
		Bottle No.:	0255

Function	Result %BAC	Time HH:MM	Temperature Simulator (°C)	Comment(s) or Error(s)
Ambient Air Blank	0.000%	16:09D		
Control 1 EC	0.101%	16:10D	34.0°C	*** TEST PASSED ***
Control 1 IR	0.101%	16:10D	34.0°C	*** TEST PASSED ***
Ambient Air Blank	0.000%	16:11D		
Control 2 EC	0.100%	16:11D	34.0°C	*** TEST PASSED ***
Control 2 IR	0.102%	16:11D	34.0°C	*** TEST PASSED ***
Ambient Air Blank	0.000%	16:12D		
Control 3 EC	0.100%	16:12D	34.0°C	*** TEST PASSED ***
Control 3 IR	0.101%	16:12D	34.0°C	*** TEST PASSED ***
Ambient Air Blank	0.000%	16:13D		

All tests within acceptable tolerance.

On this date, I installed the above indicated "NEW SOLUTION" in accordance with Alcotest 7110 operator training and procedures established by the (NJSP) Chief Forensic Scientist.

Changed By:

Last Name: JANKOWSKI

First Name: JOSEPH

MI: E

Signature:

Badge No.: 714

Date: 07/17/2010

Calibrating Unit

New Standard Solution Report

Equipment Alcotest 7110 MKIII-C Serial No.: ARTL-0005
Location: PLAINSBORO TOWNSHIP PD
Calibration File No.: 00780 Calib. Date: 06/07/2010 Calib. No.: 00014
Certification File No.: 00781 Cert. Date: 06/07/2010 Cert. No.: 00010
Linearity File No.: 00782 Lin. Date: 06/07/2010 Lin. No.: 00010
Solution File No.: 00824 Soln. Date: 08/09/2010 Soln. No.: 00105
Sequential File No.: 00824 File Date: 08/09/2010

Calibrating Unit: WET Model No.: CU-34 Serial No.: DDUF S3-0064
Control Solution %: 0.100% Expires: 12/31/2011
Solution Control Lot: 09L072 Bottle No.: 0256

Function	Result	Time	Temperature	Comment(s)
	%BAC	HH:MM	Simulator (°C)	or Error(s)
Ambient Air Blank	0.000%	08:52D		
Control 1 EC	0.100%	08:52D	34.0°C	*** TEST PASSED ***
Control 1 IR	0.101%	08:52D	34.0°C	*** TEST PASSED ***
Ambient Air Blank	0.000%	08:53D		
Control 2 EC	0.098%	08:54D	34.0°C	*** TEST PASSED ***
Control 2 IR	0.100%	08:54D	34.0°C	*** TEST PASSED ***
Ambient Air Blank	0.000%	08:55D		
Control 3 EC	0.099%	08:55D	34.0°C	*** TEST PASSED ***
Control 3 IR	0.100%	08:55D	34.0°C	*** TEST PASSED ***
Ambient Air Blank	0.000%	08:56D		

All tests within acceptable tolerance.


On this date, I installed the above indicated "NEW SOLUTION" in accordance with Alcotest 7110 operator training and procedures established by the (NJSP) Chief Forensic Scientist.

Changed By:

Last Name: TAVENER

First Name: FREDERICK

MI: W

Signature: 

Badge No.: 730

Date: 08/09/2010

Calibrating Unit

New Standard Solution Report

Equipment Alcotest 7110 MKIII-C Serial No.: ARTL-0005
Location: PLAINSBORO TOWNSHIP PD
Calibration File No.: 00780 **Calib. Date:** 06/07/2010 **Calib. No.:** 00014
Certification File No.: 00781 **Cert. Date:** 06/07/2010 **Cert. No.:** 00010
Linearity File No.: 00782 **Lin. Date:** 06/07/2010 **Lin. No.:** 00010
Solution File No.: 00830 **Soln. Date:** 08/15/2010 **Soln. No.:** 00106
Sequential File No.: 00830 **File Date:** 08/15/2010

Calibrating Unit: WET **Model No.:** CU-34 **Serial No.:** DDUF S3-0064
Control Solution %: 0.100% **Expires:** 12/31/2011
Solution Control Lot: 09L072 **Bottle No.:** 0261

Function	Result %BAC	Time HH:MM	Temperature Simulator (°C)	Comment(s) or Error(s)
Ambient Air Blank	0.000%	04:32D		
Control 1 EC	0.101%	04:32D	34.0°C	*** TEST PASSED ***
Control 1 IR	0.100%	04:32D	34.0°C	*** TEST PASSED ***
Ambient Air Blank	0.000%	04:33D		
Control 2 EC	0.100%	04:33D	34.0°C	*** TEST PASSED ***
Control 2 IR	0.102%	04:33D	34.0°C	*** TEST PASSED ***
Ambient Air Blank	0.000%	04:34D		
Control 3 EC	0.100%	04:35D	34.0°C	*** TEST PASSED ***
Control 3 IR	0.100%	04:35D	34.0°C	*** TEST PASSED ***
Ambient Air Blank	0.000%	04:35D		

All tests within acceptable tolerance.

On this date, I installed the above indicated "NEW SOLUTION" in accordance with Alcotest 7110 operator training and procedures established by the (NJSP) Chief Forensic Scientist.

Changed By:

Last Name: OLMA

First Name: BRETT

MI: J

Signature: _____

[Handwritten Signature] #35

Badge No.: 535

Date: 08/15/2010

Calibrating Unit New Standard Solution Report

Equipment	Alcotest 7110 MKIII-C	Serial No.:	ARTL-0005
Location:	PLAINSBORO TOWNSHIP PD		
Calibration File No.:	00780	Calib. Date:	06/07/2010
Certification File No.:	00781	Cert. Date:	06/07/2010
Linearity File No.:	00782	Lin. Date:	06/07/2010
Solution File No.:	00832	Soln. Date:	08/16/2010
Sequential File No.:	00832	File Date:	08/16/2010
Calibrating Unit:	WET	Model No.:	CU-34
Control Solution %:	0.100%	Serial No.:	DDUF S3-0064
Solution Control Lot:	09L072	Expires:	12/31/2011
		Bottle No.:	0257

Function	Result %BAC	Time HH:MM	Temperature Simulator (°C)	Comment(s) or Error(s)
Ambient Air Blank	0.000%	15:16D		
Control 1 EC	0.101%	15:16D	34.0°C	*** TEST PASSED ***
Control 1 IR	0.101%	15:16D	34.0°C	*** TEST PASSED ***
Ambient Air Blank	0.000%	15:17D		
Control 2 EC	0.100%	15:17D	34.0°C	*** TEST PASSED ***
Control 2 IR	0.101%	15:17D	34.0°C	*** TEST PASSED ***
Ambient Air Blank	0.000%	15:18D		
Control 3 EC	0.101%	15:19D	34.0°C	*** TEST PASSED ***
Control 3 IR	0.100%	15:19D	34.0°C	*** TEST PASSED ***
Ambient Air Blank	0.000%	15:19D		

All tests within acceptable tolerance.

On this date, I installed the above indicated "NEW SOLUTION" in accordance with Alcotest 7110 operator training and procedures established by the (NJSP) Chief Forensic Scientist.

Changed By:

Last Name: TAVENER

First Name: FREDERICK

MI: W

Signature: 

#730

Badge No.: 730

Date: 08/16/2010

Calibrating Unit

New Standard Solution Report

Equipment	Alcotest 7110 MKIII-C	Serial No.:	ARTL-0005
Location:	PLAINSBORO TOWNSHIP PD		
Calibration File No.:	00780	Calib. Date:	06/07/2010
Certification File No.:	00781	Cert. Date:	06/07/2010
Linearity File No.:	00782	Lin. Date:	06/07/2010
Solution File No.:	00848	Soln. Date:	09/10/2010
Sequential File No.:	00848	File Date:	09/10/2010
Calibrating Unit:	WET	Model No.:	CU-34
Control Solution %:	0.100%	Serial No.:	DDUF S3-0064
Solution Control Lot:	09L072	Expires:	12/31/2011
		Bottle No.:	0258

Function	Result %BAC	Time HH:MM	Temperature Simulator (°C)	Comment(s) or Error(s)
Ambient Air Blank	0.000%	19:04D		
Control 1 EC	0.100%	19:04D	34.0°C	*** TEST PASSED ***
Control 1 IR	0.101%	19:04D	34.0°C	*** TEST PASSED ***
Ambient Air Blank	0.000%	19:05D		
Control 2 EC	0.099%	19:05D	34.0°C	*** TEST PASSED ***
Control 2 IR	0.101%	19:05D	34.0°C	*** TEST PASSED ***
Ambient Air Blank	0.000%	19:06D		
Control 3 EC	0.100%	19:07D	34.0°C	*** TEST PASSED ***
Control 3 IR	0.102%	19:07D	34.0°C	*** TEST PASSED ***
Ambient Air Blank	0.000%	19:07D		

All tests within acceptable tolerance.

On this date, I installed the above indicated "NEW SOLUTION" in accordance with Alcotest 7110 operator training and procedures established by the (NJSP) Chief Forensic Scientist.

Changed By:

Last Name: TAVENER

First Name: FREDERICK

MI: W

Signature: SWT-Flan

Badge No.: 730

Date: 09/10/2010

Calibrating Unit

New Standard Solution Report

Equipment Alcotest 7110 MKIII-C Serial No.: ARTL-0005
Location: PLAINSBORO TOWNSHIP PD
Calibration File No.: 00780 Calib. Date: 06/07/2010 Calib. No.: 00014
Certification File No.: 00781 Cert. Date: 06/07/2010 Cert. No.: 00010
Linearity File No.: 00782 Lin. Date: 06/07/2010 Lin. No.: 00010
Solution File No.: 00856 Soln. Date: 10/13/2010 Soln. No.: 00109
Sequential File No.: 00856 File Date: 10/13/2010

Calibrating Unit: WET Model No.: CU-34 Serial No.: DDUF S3-0064
Control Solution %: 0.100% Expires: 12/31/2011
Solution Control Lot: 09L072 Bottle No.: 0259

Function	Result %BAC	Time HH:MM	Temperature Simulator (°C)	Comment(s) or Error(s)
Ambient Air Blank	0.000%	16:21D		
Control 1 EC	0.101%	16:21D	34.0°C	*** TEST PASSED ***
Control 1 IR	0.101%	16:21D	34.0°C	*** TEST PASSED ***
Ambient Air Blank	0.000%	16:22D		
Control 2 EC	0.100%	16:22D	34.0°C	*** TEST PASSED ***
Control 2 IR	0.100%	16:22D	34.0°C	*** TEST PASSED ***
Ambient Air Blank	0.000%	16:23D		
Control 3 EC	0.100%	16:24D	34.0°C	*** TEST PASSED ***
Control 3 IR	0.100%	16:24D	34.0°C	*** TEST PASSED ***
Ambient Air Blank	0.000%	16:24D		

All tests within acceptable tolerance.

On this date, I installed the above indicated "NEW SOLUTION" in accordance with Alcotest 7110 operator training and procedures established by the (NJSP) Chief Forensic Scientist.

Changed By:

Last Name: TAVENER

First Name: FREDERICK

MI: W

Signature: SGT F Taver

#730

Badge No.: 730

Date: 10/13/2010

Calibrating Unit

New Standard Solution Report

Equipment Alcotest 7110 MKIII-C Serial No.: ARTL-0005
Location: PLAINSBORO TOWNSHIP PD
Calibration File No.: 00780 Calib. Date: 06/07/2010 Calib. No.: 00014
Certification File No.: 00781 Cert. Date: 06/07/2010 Cert. No.: 00010
Linearity File No.: 00782 Lin. Date: 06/07/2010 Lin. No.: 00010
Solution File No.: 00870 Soln. Date: 11/12/2010 Soln. No.: 00110
Sequential File No.: 00870 File Date: 11/12/2010

Calibrating Unit: WET Model No.: CU-34 Serial No.: DDUF S3-0064
Control Solution %: 0.100% Expires: 12/31/2011
Solution Control Lot: 09L072 Bottle No.: 0262

Function	Result %BAC	Time HH:MM	Temperature Simulator (°C)	Comment(s) or Error(s)
Ambient Air Blank	0.000%	12:18S		
Control 1 EC	0.102%	12:19S	34.0°C	*** TEST PASSED ***
Control 1 IR	0.101%	12:19S	34.0°C	*** TEST PASSED ***
Ambient Air Blank	0.000%	12:19S		
Control 2 EC	0.100%	12:20S	34.0°C	*** TEST PASSED ***
Control 2 IR	0.101%	12:20S	34.0°C	*** TEST PASSED ***
Ambient Air Blank	0.000%	12:21S		
Control 3 EC	0.100%	12:21S	34.0°C	*** TEST PASSED ***
Control 3 IR	0.101%	12:21S	34.0°C	*** TEST PASSED ***
Ambient Air Blank	0.000%	12:22S		

All tests within acceptable tolerance.

On this date, I installed the above indicated "NEW SOLUTION" in accordance with Alcotest 7110 operator training and procedures established by the (NJSP) Chief Forensic Scientist.

Changed By:

Last Name: DELL BENI

First Name: JOSEPH

MI: M

Signature: RI Joseph Beni

543

Badge No.: 543

Date: 11/12/2010

Calibrating Unit

New Standard Solution Report

Equipment

Alcotest 7110 MKIII-C
Location: PLAINSBORO TOWNSHIP PD
Calibration File No.: 00879
Certification File No.: 00880
Linearity File No.: 00881
Solution File No.: 00882
Sequential File No.: 00882

Calib. Date: 11/22/2010
Cert. Date: 11/22/2010
Lin. Date: 11/22/2010
Soln. Date: 11/22/2010
File Date: 11/22/2010

Serial No.: ARTL-0005
Calib. No.: 00015
Cert. No.: 00011
Lin. No.: 00011
Soln. No.: 00111

Calibrating Unit: WET
Control Solution %: 0.100%
Solution Control Lot: 09L072

Model No.: CU-34

Serial No.: DDUF S3-0064
Expires: 12/31/2011
Bottle No.: 0268

Function	Result	Time	Temperature	Comment(s)
	%BAC	HH:MM	Simulator (°C)	or Error(s)
Ambient Air Blank	0.000%	15:13S		
Control 1 EC	0.100%	15:13S	34.0°C	*** TEST PASSED ***
Control 1 IR	0.098%	15:13S	34.0°C	*** TEST PASSED ***
Ambient Air Blank	0.000%	15:14S		
Control 2 EC	0.099%	15:14S	34.0°C	*** TEST PASSED ***
Control 2 IR	0.098%	15:14S	34.0°C	*** TEST PASSED ***
Ambient Air Blank	0.000%	15:15S		
Control 3 EC	0.098%	15:16S	34.0°C	*** TEST PASSED ***
Control 3 IR	0.097%	15:16S	34.0°C	*** TEST PASSED ***
Ambient Air Blank	0.000%	15:16S		

All tests within acceptable tolerance.

On this date, I installed the above indicated "NEW SOLUTION" in accordance with Alcotest 7110 operator training and procedures established by the (NJSP) Chief Forensic Scientist.

TEMPERATURE PROBE SERIAL #DDSF P2-0396

ND

Changed By:

Last Name: DENNIS

First Name: MARC

MI: W.

Signature:

TAL I [Signature] #5925

Badge No.: 5925

Date: 11/22/2010

Calibrating Unit

New Standard Solution Report

Equipment Alcotest 7110 MKIII-C Serial No.: ARTL-0005
Location: PLAINSBORO TOWNSHIP PD
Calibration File No.: 00879 Calib. Date: 11/22/2010 Calib. No.: 00015
Certification File No.: 00880 Cert. Date: 11/22/2010 Cert. No.: 00011
Linearity File No.: 00881 Lin. Date: 11/22/2010 Lin. No.: 00011
Solution File No.: 00891 Soln. Date: 12/11/2010 Soln. No.: 00112
Sequential File No.: 00891 File Date: 12/11/2010

Calibrating Unit: WET Model No.: CU-34 Serial No.: DDUF S3-0064
Control Solution %: 0.100% Expires: 12/31/2011
Solution Control Lot: 09L072 Bottle No.: 0263

Function	Result %BAC	Time HH:MM	Temperature Simulator (°C)	Comment(s) or Error(s)
Ambient Air Blank	0.000%	08:15S		
Control 1 EC	0.099%	08:15S	34.0°C	*** TEST PASSED ***
Control 1 IR	0.098%	08:15S	34.0°C	*** TEST PASSED ***
Ambient Air Blank	0.000%	08:16S		
Control 2 EC	0.098%	08:17S	34.0°C	*** TEST PASSED ***
Control 2 IR	0.098%	08:17S	34.0°C	*** TEST PASSED ***
Ambient Air Blank	0.000%	08:18S		
Control 3 EC	0.097%	08:18S	34.0°C	*** TEST PASSED ***
Control 3 IR	0.099%	08:18S	34.0°C	*** TEST PASSED ***
Ambient Air Blank	0.000%	08:19S		

All tests within acceptable tolerance.

On this date, I installed the above indicated "NEW SOLUTION" in accordance with Alcotest 7110 operator training and procedures established by the (NJSP) Chief Forensic Scientist.

Changed By:

Last Name: DELL BENI

First Name: JOSEPH

MI: M

Signature: P. J. Joseph 8/11/10

543

Badge No.: 543

Date: 12/11/2010

Calibrating Unit

New Standard Solution Report

Equipment Alcotest 7110 MKIII-C Serial No.: ARTL-0005
Location: PLAINSBORO TOWNSHIP PD
Calibration File No.: 00879 **Calib. Date:** 11/22/2010 **Calib. No.:** 00015
Certification File No.: 00880 **Cert. Date:** 11/22/2010 **Cert. No.:** 00011
Linearity File No.: 00881 **Lin. Date:** 11/22/2010 **Lin. No.:** 00011
Solution File No.: 00902 **Soln. Date:** 12/30/2010 **Soln. No.:** 00113
Sequential File No.: 00902 **File Date:** 12/30/2010

Calibrating Unit: WET **Model No.:** CU-34 **Serial No.:** DDUF S3-0064
Control Solution %: 0.100% **Expires:** 12/31/2011
Solution Control Lot: 09L072 **Bottle No.:** 0264

Function	Result %BAC	Time HH:MM	Temperature Simulator (°C)	Comment(s) or Error(s)
Ambient Air Blank	0.000%	14:04S		
Control 1 EC	0.100%	14:04S	34.0°C	*** TEST PASSED ***
Control 1 IR	0.098%	14:04S	34.0°C	*** TEST PASSED ***
Ambient Air Blank	0.000%	14:05S		
Control 2 EC	0.098%	14:06S	34.0°C	*** TEST PASSED ***
Control 2 IR	0.099%	14:06S	34.0°C	*** TEST PASSED ***
Ambient Air Blank	0.000%	14:07S		
Control 3 EC	0.098%	14:07S	34.0°C	*** TEST PASSED ***
Control 3 IR	0.098%	14:07S	34.0°C	*** TEST PASSED ***
Ambient Air Blank	0.000%	14:08S		

All tests within acceptable tolerance.


On this date, I installed the above indicated "NEW SOLUTION" in accordance with Alcotest 7110 operator training and procedures established by the (NJSP) Chief Forensic Scientist.

Changed By:

Last Name: DELL BENI

First Name: JOSEPH

MI: M

Signature: 

Badge No.: 543

Date: 12/30/2010